

(19)



Europäisches Patentamt
European Patent Office
Office européen des brevets



(11) Publication number:

0 422 872 A2

(12)

EUROPEAN PATENT APPLICATION

(21) Application number: 90311003.9

(51) Int. Cl.⁵: **C12Q 1/68, C07H 21/04**

(22) Date of filing: 08.10.90


(30) Priority: 12.10.89 US 420577

(43) Date of publication of application:
17.04.91 Bulletin 91/16(84) Designated Contracting States:
AT BE CH DE DK ES FR GB GR IT LI LU NL SE(71) Applicant: **GENE-TRAK SYSTEMS**
31 New York Avenue
Framingham Massachusetts 01701(US)(72) Inventor: **Weisburg, William G.**
3 Jillson Circle
Milford, Massachusetts 01757(US)

Inventor: **Sogin, Mitchell I.**
1371 GlenCoe Street
Denver, Colorado 80220(US)
Inventor: **Barns, Susan M.**
3 Maple Street
Hopkinton, Massachusetts 01748(US)
Inventor: **Pelletier, Dale A.**
58 Walnut Street, Apartment 3
Somerville, Massachusetts 02143(US)

(74) Representative: **Sheard, Andrew Gregory et al**
Kilburn & Strode 30, John Street
London WC1N 2DD(GB)(54) **Nucleic acid probes and methods for detecting fungi.**

(57) Nucleic acid probes are described for detecting fungi capable of causing fungal septicemia or capable of causing food spoilage. The preferred probes are complementary to ribonucleic acid sequences found in numerous fungi and absent in animal or plant genomes. As such, these probes can detect the rRNA, rDNA, or polymerase chain reaction amplification products from the majority of fungal species. The detection of etiological agents of human fungemia, the clinical diagnosis of this disease and the direct evaluation of food or beverage fungal content utilizing rRNA or rDNA probes is now possible.

FIGURE

NUCLEIC ACID PROBES AND METHODS FOR DETECTING FUNGI

Field of the Invention

This invention relates to the detection of fungal organisms. More specifically, it provides nucleic acid probes and compositions along with methods for their use for the specific detection of yeasts and molds in clinical, food, environmental and other samples.

The fungi are a diverse collection of cell-wall enclosed eukaryotes either saprophytic or parasitic and may be morphologically described as yeasts, molds, mushrooms, or by other names. They are ubiquitous organisms, mostly innocuous, sometimes used for commercial purposes, and occasionally pathogenic.

The pathogenic fungi are included within the domain of medical mycology. This medical field recognizes categories of fungal pathogens (see Rippon, J.W., Medical Mycology, Saunders Co., Philadelphia, 1988, for example) including superficial, cutaneous, subcutaneous, and systemic infection. By far the most serious pathology caused by the fungi that clinicians face are the systemic infections. Deep tissue and systemic fungemia claim high mortality rates, particularly among immune comprised populations.

Among the fungi capable of causing systemic fungemia, there is a dichotomy between the so-called "pathogenic" fungi and the "opportunistic" fungi. It is a deceptive nomenclature; the opportunists are the killers, and the pathogenic fungi are often self-limiting. The pathogenic fungi include *Coccidioides immitis*, *Histoplasma capsulatum*, *Blastomyces dermatitidis*, *Paracoccidioides brasiliensis*, and the subcutaneous pathogen, *Sporothrix schenckii*. The important opportunistic fungi include the *Candida* -- particularly *C. albicans*, *C. tropicalis*, *C. parapsilosis*, and *Torulopsis (Candida) glabrata* -- *Cryptococcus neoformans*, members of the genus *Aspergillus*, and to a lesser extent, practically any fungus that can survive at host physiological temperatures.

Clinical diagnosis and treatment of systemic fungemia suffers several shortcomings compared to bacterial septicemia (which of ten occurs in the same immune-deficient population). First, antifungal chemotherapy is more toxic to the patient than analogous antibacterial chemotherapy. As a result, clinicians desire a more reliable demonstration of fungemia before prescribing antifungal agents. Second, fungemic patients have a poor prognosis, unless diagnosed early in infection. Third, fungi generally grow slower than the major bacteremic organisms, and consequently diagnosis requiring an

in vitro culture step is time consuming. And fourth, some of the fungi (again in diagnoses requiring in vitro cultivation) will not yield colonies on synthetic media for weeks, if at all. All of these factors, plus the fact that a wide array of fungi are potential systemic pathogens, point to the need for a direct method of fungal detection inclusive for virtually all fungi.

It is an aspect of the present invention to provide nucleic acid probes capable of detecting fungi.

It is another aspect of the present invention to provide nucleic acid probes which can hybridize to target regions which can be rendered accessible to probes under normal assay conditions.

It is yet another aspect to provide nucleic acid probes to fungal rRNA sequences useful as the basis for rapid diagnostic assays for assessing the presence of these organisms in a clinical sample.

While Kohne et al. (Biophysical Journal 8:1104-1118, 1968) discuss one method for preparing probes to rRNA sequences, they do not provide the teaching necessary to make probes to detect fungi. Pace and Campbell (Journal of Bacteriology 107:543-547, 1971) discuss the homology of ribosomal ribonucleic acids from diverse bacterial species and a hybridization method for quantitating such homology levels. Similarly, Sogin, Sogin and Woese (Journal of Molecular Evolution 1:173-184, 1972) discuss the theoretical and practical aspects of using primary structural characterization of different ribosomal RNA molecules for evaluating phylogenetic relationships. Pox, Pechman and Woese (International Journal of Systematic Bacteriology 27:44-57, 1977) discuss the comparative cataloging of 16S ribosomal RNAs as an approach to prokaryotic systematics. These references, however, fail to relieve the deficiency of Kohne's teaching with respect to fungi, and in particular, do not provide specific probes useful in assays for detecting fungemia or its etiological agents, a broad spectrum of yeast and molds.

Hogan, et al (International Patent Application, Publication Number WO 88/03957) describe four putative fungal specific probes. None of them appear widely inclusive for fungi, nor are they related to the probes of the present invention.

Ribosomes are of profound importance to all organisms because they serve as the only known means of translating genetic information into cellular proteins, the main structural and catalytic elements of life. A clear manifestation of this importance is the observation that all cells have ribosomes.

Bacterial ribosomes contain three distinct RNA

molecules which, at least in *Escherichia coli*, are referred to as 5S, 16S and 23S rRNAs. In eukaryotic organisms, there are four distinct rRNA species, generally referred to as 5S, 18S, 28S, and 5.8S. These names historically are related to the size of the RNA molecules, as determined by their sedimentation rate. In actuality, however, ribosomal RNA molecules vary substantially in size between organisms. Nonetheless, 5S, 18S, 28S, and 5.8S rRNA are commonly used as generic names for the homologous RNA molecules in any eukaryote, and this convention will be continued herein.

It is another aspect of the present invention to provide nucleic acid probes complimentary to unique nucleic acid sequences within the 18S ribosomal ribonucleic acid (rRNA) of fungal pathogens.

As used herein, probe(s) refer to synthetic or biologically produced nucleic acids (DNA or RNA) which, by design or selection, contain specific nucleotide sequences that allow them to hybridize under defined predetermined stringencies, specifically (i.e., preferentially, see next paragraph) to target nucleic acid sequences. In addition to their hybridization properties, probes also may contain certain constituents that pertain to their proper or optimal functioning under particular assay conditions. For example, probes may be modified to improve their resistance to nuclease degradation (e.g. by end capping), to carry detection ligands (e.g. fluorescein, 32-P, biotin, etc.), or to facilitate their capture onto a solid support (e. g., polydeoxyadenosine "tails"). Such modifications are elaborations on the basic probe function which is its ability to usefully discriminate between target and non-target organisms in a hybridization assay.

Hybridization traditionally is understood as the process by which, under predetermined reaction conditions, two partially or completely complementary strands of nucleic acid are allowed to come together in an antiparallel fashion (one oriented 5' to 3', the other 3' to 5') to form a double-stranded nucleic acid with specific and stable hydrogen bonds, following explicit rules pertaining to which nucleic acid bases may pair with one another. The high specificity of probes relies on the low statistical probability of unique sequences occurring at random as dictated by the multiplicative product of their individual probabilities. These concepts are well understood by those skilled in the art.

The stringency of a particular set of hybridization conditions is determined by the base composition of the probe/target duplex, as well as by the level and geometry of mispairing between the two nucleic acids.

Stringency may also be governed by such reaction parameters as the concentration and type of ionic species present in the hybridization solu-

tion, the types and concentrations of denaturing agents present, and the temperature of hybridization. Generally, as hybridization conditions become more stringent, longer probes are preferred if stable hybrids are to be formed. As a corollary, the stringency of the conditions under which a hybridization is to take place (e. g., based on the type of assay to be performed) will dictate certain characteristics of the preferred probes to be employed. Such relationships are well understood and can be readily manipulated by those skilled in the art.

As a general matter, dependent upon probe length, such persons understand stringent conditions to mean approximately 35°C-65°C in a salt solution of approximately 0.9 molar NaCl.

All references made herein are fully incorporated by reference.

In accordance with the various principles and aspects of the present invention, there are provided nucleic acid probes and probe sets comprising deoxyribonucleic acid (DNA) or ribonucleic acid (RNA) sequences which hybridize, under specific conditions, to the ribosomal RNA molecules (rRNA), specifically 18S rRNA molecules, or rRNA genes (rDNA) of fungi but which do not hybridize, under the same conditions, to the rRNA or rDNA of bacteria or the host or environmental matrix which may be present in test samples. The probes of the present invention now permit the development of a valuable nucleic acid hybridization assay for the specific detection of fungemia or its etiological agents. This assay may advantageously be used to test for yeasts and molds in clinical samples of blood, urine, cerebrospinal fluid, skin biopsy, saliva, synovial fluid, sputum, bronchial wash, bronchial lavage, or other tissue or fluid samples from human patients or veterinary subjects.

The probes of the present invention also provide the basis for the development of valuable nucleic acid hybridization assays capable of detecting yeasts and molds associated with food spoilage. Most preferred probes of the present invention can hybridize to a diverse collection of fungi while not cross-reacting, at predetermined conditions, with meats, dairy products, grains, nuts, juices, and other commercial food matrices.

Nucleic acid hybridization based assays have been discovered to impart enhanced performance capabilities with respect to most currently available, microbiological or immunological methods for detection of fungi in test samples, generally including:

- a) increased sensitivity; i.e., the ability to detect yeast or mold in a given sample more frequently;
- b) potentially significant reductions in assay cost due to the use of inexpensive reagents and reduced labor;
- c) accurate identification of even biochemically

unusual strains of the target organism, or isolates with dramatically different antigenic properties;

d) direct assay for the presence of the yeast or mold and consequent potential to quantify the etiological agents;

e) direct testing allows the monitoring of the efficacy of an antifungal regime; and f) potentially significant reductions in the exposure of laboratory technologists to bodily fluid specimens harboring infectious agents.

It has been discovered that other advantages incurred by directing the probes of the present invention against rRNA include the fact that the rRNAs detected constitute a significant component of cellular mass. Although estimates of cellular ribosome content vary, actively growing fungal cells may contain upwards of 100,000 ribosomes per cell, and therefore 100,000 copies of each of the rRNAs (present in a 1:1:1:1 stoichiometry in ribosomes). In contrast, other potential cellular target molecules such as genes or RNA transcripts thereof, are less ideal since they are present in much lower abundance. A further unexpected advantage is that the rRNAs (and the genes specifying them) appear not to be subject to lateral transfer between contemporary organisms. Thus, the rRNA primary structure provides an organism-specific molecular target, rather than a gene-specific target as would likely be the case, for example of a plasmid-borne gene or product thereof which may be subject to lateral transmission between contemporary organisms.

The discovery that probes could be generated with the extraordinary inclusivity and exclusivity characteristics of those of the present invention with respect to the detection of virtually all fungal organism, without necessarily incurring cross-reactivity to animal, plant, or bacterial genomes was unpredictable and unexpected.

Further understanding of the principles and aspects of the present invention may be made by reference to the tables wherein:

Tables 1, 2, and 3 display the hybridization behavior of fifteen probes toward a panel of clinically and environmentally representative fungal species. Additional fungi were added to the panel in order to represent the breadth of known fungal taxa. Approximately eighty fungal species are represented, and the highest prevalence pathogens are represented by numerous strains. In addition, nucleic acids from a variety of non-fungal organisms are included for comparison including RNAs from human, wheat, normal human stool, and two ubiquitous bacterial species. Those skilled in the art understand that bacteria are so evolutionarily distant as to not generally cross-react with the types of probes described herein. It will be further

recognized that the sequence variation among vertebrate animals and among higher plants is sufficiently narrow that an individual sample, such as wheat, has high predictive value.

All species on the panel are represented by 100 ng of purified, denatured RNA. Probes were 32-Phosphorous labelled, hybridized to panels under standard conditions, at the temperatures indicated, and autoradiographically evaluated. "+" represents strong hybridization signal after three hours exposure, "+-" is a weak signal, "+--" is virtually absent, and "-" is indicative of no hybridization of probe to target. "NT" indicates that a particular probe was not tested against the designated strain.

Still further understanding may be had by study of the accompanying Figure which shows a schematic representation of a dual probe capture/detector assay.

Probe Development Strategy

The first step taken in the development of the probes of the present invention involved identification of regions of the 18S rRNA which potentially could serve as target sites for fungal specific nucleic acid probes. This entailed finding sites which are:

- 1) highly conserved (few nucleotide changes, deletions, or insertions) among the fungal rRNA sequences, and
- 2) substantially different in non-fungal (bacteria, human, or plant) rRNA sequences.

For this analysis, precise alignments of available 18S rRNA sequences were developed. A number of 18S rRNA sequences were determined as part of this effort. Such nucleotide sequences were determined by standard laboratory protocols either by cloning and sequencing of genes specifying rRNAs or by direct sequencing of the rRNAs themselves using reverse transcriptase (Lane, et al, 1985, Proceedings of the National Academy of Sciences, USA 82:6955-6959).

A computer algorithm, operating on the aligned set of 18S rRNA sequences, was used to identify regions of greatest similarity among the fungi. Nucleic acid probes to such regions will hybridize most widely among diverse fungi. Additional information was gained by comparing these fungal conserved regions to known 18S rRNA sequences from human, rat, mouse, corn, soy, rice, bacteria, protozoa, algae, etc.

Fifteen probes were identified, based on these analyses. The discovery of a specific type of non-fungal cross-reactivity does not necessarily render a probe uninteresting. For example, cross-reaction of a probe with the plant kingdom does not detract

from the probe's usefulness in screening vertebrate blood for fungemia. Other probes described herein are known to cross-react, but they were designed to be employed in dual probe assays (see Figure and Example 2).

Probes may be synthesised by any suitable chemical or biological method. In particular, DNA probes may be synthesised by automated phosphoramidite chemistry using cyanoethyl phosphoramidites (Beaucage and Carruthers *Tetrahedron Letters* 24, 245 (1981)). RNA probes can be made by transcription of the corresponding DNA sequences. screening vertebrate blood for fungemia. Other probes described herein are known to cross-react, but they were designed to be employed in dual probe assays (see the Figure and Example 2).

Description of the Probes

As indicated, the above probe selection strategy yielded fifteen probes useful for hybridization to fungi in samples comprising:

PROBE 1417: 5' - TGTCTGGACCTGGTGAGTTTCCCCGTG-3' (SEQ ID1)
 PROBE 1418: 5' - TGTCTGGACCTGGTGAGTTTCCCCGTGTTGAGTCAAATT-3' (SEQ ID2)
 PROBE 1415: 5' - TCCTGGTTAAGGGATTAAATTGTACTCATTCCAATT-3' (SEQ ID3)
 PROBE 1416: 5' - TCCTCGTTAAGGTATTTACATTGTACTCATTCCAATT-3' (SEQ ID4)
 PROBE 1G707: 5' - TCCTCGTTAAGGTOTTTAAATTGTACTCATTCCAATT-3' (SEQ ID5)
 PROBE 1542: 5' - AACTAAGAACGGCCATGCACCACCAT-3' (SEQ ID6)
 PROBE 1545: 5' - TGGTGCCCTTCGGTCAATTTCTTTAATCAGCCTTGCG-3' (SEQ ID7)
 PROBE 1814: 5' - TCGCTGGCGCAAGGCCATGCGATTGAGAGGTTATTATGAATCATCAG-3' (SEQ ID8)
 PROBE 1816: 5' - CAAGCTGATGACTTGTGCTTACTAGGGATT-3' (SEQ ID9)
 PROBE 1857: 5' - TCGGCATAGTTTGTGGTTAAGACTACGACGGTATCTT-3' (SEQ ID10)
 PROBE 1813: 5' - AAATGCTTTCGCAGTAGTTGGTCTT-3' (SEQ ID11)
 PROBE 1860: 5' - AAATGGCGCAGTAGTTGGTCT-

TCGGTAAATCCAAGAATTTACCTT-3' (SEQ ID12)
 PROBE 1812: 5' - ACGTCCTATTTTATTATTCCATGCTAAT-3' (SEQ ID 13)
 PROBE 1858: 5' - AAGTCATATTTTCATTATTCCATGCTAACT-3' (SEQ ID14)
 PROBE 1859: 5' - TCGTCGAGTTATGTTATTCCATGCAAAT-3' (SEQ ID15)

The specific behaviors of the aforementioned probes are dependent to a significant extent on the assay format in which they are employed. Conversely, the assay format will dictate certain of the optimal features of particular probes. The "essence" of the probes of the invention is not to be construed as restricted to the specific string of nucleotides in the named probes. For example, the length of these particular oligonucleotides was optimized for use in the dot blot assay (and certain other anticipated assays) described below. It is well known to those skilled in the art that optimal probe length will be a function of the stringency of the hybridization conditions chosen and hence the length of the instant probes may be altered accordingly. Also, in considering sets comprised of more than one probe, it is desirable that all probes behave in a compatible manner in any particular format in which they are employed. Thus, the exact length of a particular probe will to a certain extent, reflect its specific intended use.

The probes of the present invention are useful as oligonucleotide probes and can also be incorporated into larger polynucleotides of either ribonucleic acid or deoxyribonucleic acid. Sequences complimentary to the probes described herein can be used as probes to rRNA genes. The preferred probes or their compliments also can be employed as chain elongation initiators for polymerase chain reaction, sequencing or other applications.

Two additional preferred probes useful in this regard comprise:

Probe/Primer 936: 5' -(CCGAATTCGTCGACAAC)-CTGGTTGATCCTGCCAGT-3' (SEQ ID 16)

Probe/Primer 935: 5' -(CCCGGGATCGAAGCT)-TGATCCTTCTGCAGGTTACCTAC-3' (SEQ ID 17)

Probe/Primer 936 is designed to hybridize to the 18S rDNA gene strand complimentary to fungal 18S rRNA. Oligonucleotides 935 and 936 are designed for and most preferred for use in assays employing amplification, by the polymerase chain reaction method, of almost the entire 18S rRNA gene (rDNA) of fungi and relatives. The target specific "essence" of these two probe/primers resides in the portions of these oligonucleotides not included within the parenthesis. The nucleotides within the parentheses are preferably included since they add useful restriction endonuclease rec-

ognition (cloning) sites to the amplified products.

Probe Behavior During Hybridization

The experimental specificity of the preferred probes, as further documented in Example 1 and Tables 1, 2, and 3, may be summarized as follows: Probe 1417: 100% inclusive for tested fungi, with negligible cross-reactivity to human RNA at 60°C. At 65°C (hybridization temperature) signal is decreased for 2 of the 171 fungi, all others still strongly hybridize, and the human cross-reactivity is removed. Strong hybridization to wheat rRNA is evident.

Probe 1418: 100% inclusive for all tested fungi, with no cross-reactivity to human RNA.

Probe 1417: Is a subsequence of 1418, that is, it is a shorter version of the same probe.

Probe 1415: Inclusive for a subset of fungi, not including any *Candida* yeasts, but including the important pathogenic yeast, *Cryptococcus*.

Probe 1416: Inclusive of all of the tested strains from the genus *Candida* (*Torulopsis*) except for the species *Yarrowia* (*Candida*) *lipolytica*. Also inclusive for *Hansenula*, *Metschnikowia*, and *Saccharomyces* --all close evolutionary relatives of *Candida* yeasts. One *Penicillium* species yields a weak signal with this probe in this assay format.

Probe IG707: 100% inclusive for *Yarrowia lipolytica*. Combined with probe 1416, these two are fully inclusive for the *Candidas* tested.

Probes 1415, 1416, and IG707: are homologous, that is, they all hybridize to the same region of the 18S rRNA.

Probe 1542: 100% inclusive for all tested fungi plus human RNA. Designed as a companion probe for 1417 or 1418 in a dual probe (sandwich type) hybridization scheme.

Probe 1545: 100% inclusive for all tested fungi plus human RNA. Designed as a companion probe for 1417 or 1418 in a dual probe (sandwich type) hybridization scheme.

Probe 1814: Broadly fungal inclusive, particularly under 50°C hybridization conditions.

Probe 1816: Broadly fungal inclusive at 50°C.

Probe 1857: Broadly inclusive with slight hybridization to non-fungi at 50°C.

Probe 1813: Broadly inclusive at 50°C.

Probe 1860: Broadly inclusive at 50°C or 60°C hybridization.

Probe 1813: Is a subsequence of Probe 1860.

Probe 1812: Very broadly inclusive at 50°C, with no cross-reactivity to human or wheat germ RNA.

Probe 1858: Designed to hybridize to *Zygomycetes*, and thus complements the hybridizing behavior of Probe 1812.

Probe 1859: Designed to hybridize to *Yarrowia*

lipolytica, and thus complements the hybridizing behavior of Probe 1812.

Probes 1812, 1858, and 1859: Are a homologous set--that is they all hybridize to an identical location on the 18S rRNA molecule. As a set, they fail to hybridize strongly to only two strains (see Tables 2 and 3).

Non-homologous probes, such as Probes 1857 and 1860 are designed to be used together in dual probe assays as described in example 2.

Probe/primers 935 and 936 have been used to amplify 18S rDNA from all fungal taxa tested, including *Aspergillus*, *Candida*, *Penicillium*, *Cryptococcus*, and *Blastomyces*.

Example 1 Dot-Blot Analysis of Probe Hybridization Behavior

Dot-blot analysis, in accordance with well known procedures, involves immobilizing a nucleic acid or a population of nucleic acids on a filter such as nitrocellulose, nylon, or other derivatized membranes which can readily be obtained commercially, specifically for this purpose. Either DNA or RNA can be easily immobilized on such a filter and subsequently can be probed or tested for hybridization under any of a variety of conditions (i.e., stringencies) with nucleotide sequences or probes of interest. Under stringent conditions, probes whose nucleotide sequences have greater complementarity to the target will exhibit a higher level of hybridization than probes containing less complementarity.

Probes of the present invention were tested in a dot-blot format. One hundred nanograms of target RNA, purified by phenol extraction and centrifugation through cesium trifluoroacetate gradients, was denatured and spotted on a nylon membrane. Probes were isotopically labelled with the addition of a 32-Phosphorous moiety to the 5' end of the oligonucleotide. Hybridization of probes occurred, at temperatures indicated, in the presence of 1.08 M sodium chloride, 60 mM sodium phosphate, and 6 mM ethylenediamine tetraacetic acid, pH 7.4. Unhybridized probe was removed by washing at a salt concentration one-third of the hybridization condition. The filters were exposed to X-ray film and the intensity of hybridization signals was evaluated after three hours of exposure.

Tables 1, 2, and 3 summarize the behavior of the probes as tested by the above procedure and documents the specificity summarized above.

Example 2: Dual Probe Hybridization

In actual practice, many applications of these

probes would employ a pair of probes being used simultaneously in a "sandwich" hybridization scheme of "capture" probe and "detector" probe as shown in figure. The capture probe¹² ideally would be a bifunctional polynucleotide manufactured by adding a homopolymeric 3' tail to a probe with high target specificity. The tail would, in turn, hybridize to the complementary homopolymer¹¹ on a solid surface¹⁰, such as a glass bead or a filter disc. Hybridization of the capture probe¹² to its target¹⁵, in this case Fungal spirochete 18S rRNA, would complex the target¹⁵ with the solid support¹⁰. The detector probe¹³, advantageously also with some degree of specificity, would be part of a preferred detection scheme relying on radioactivity, fluorescence, chemiluminescence, color, etc. (detection moiety¹⁴) which would report the presence of the entire hybridization complex.

Example 3 Clinical Diagnosis of Fungal Septicemia from Blood, Sputum, or Cerebrospinal Fluid Sample

The clinical sample is ideally processed so as to liberate the total nucleic acid content such as by sonication, vortexing with glass beads, detergent lysis using an agent such as SDS or by chemical treatment. Alternatively, fungal cells may be partially purified by, for example, the Dupont Isolator System, followed by cell lysis. The sample, containing disrupted fungi is then incubated in the presence of capture probe, detector probe, and ideally magnetic particle beads which have been derivatized with oligo-Thymidine (see also Example 2) in a chaotropic buffer such as guanidinium isothiocyanate described by Gillespie et al USSN 299,150.

If yeast or mold 18S rRNA target molecules are present, a Bead + Capture Probe + Target + Detector Probe hybridization complex is formed. The exterior presence of a magnet near the bottom of the reaction tube will cause the magnetic particle - hybridization complex to adhere to the interior side of the tube thereby advantageously enabling removal of the unreacted components such as sample matrix, unbound probe, etc. Repeated rehydration and denaturation of the bead-probe-target complex would enable significant background reduction (as more fully described in Collins et al, USSN 922,155, EPA 87309308.2 and USSN 136,920, EPA 88312135.2). In this example, final detection could entail spotting the beads on membrane and assaying by autoradiography.

For such assays, the following capture and detector probe combinations are examples of the preferred pairs:

Probes 1417 + 1542, Probes 1417 + 1545, Probes

1418 + 1542, Probes 1418 + 1545, Probes
1416 + 1812, Probes 1812 + 1860, Probes
1857 + 1860.

Example 4 Clinical Diagnosis of Fungal Infection from Human Sample Employing Polymerase Chain Reaction Amplification of Fungal rDNA

Sample processing such as provided in Example 3 is ideally designed so as to yield DNA. The DNA is further treated to make it single stranded (e.g. by melting) in preparation for polymerase chain reaction ("PCR") amplification. Probe/Primer 936 and Probe/Primer 935 are ideally employed in conjunction with the clinical sample in the standard PCR procedures. Resultant material may then be suitably assayed utilizing the "sandwich" hybridization procedures of Example 2 with any of the probes described herein. The polymerase chain reaction can, itself, be made highly specific by employing Probe/Primer 936 in conjunction with, for example, Probe 1812. Detection is advantageously accomplished using Probe 1814 for capture and Probes 1415 and 1416 for detection.

Example 5 In situ Hybridization as a Cytological Stain

The probes of the present invention can also be advantageously employed as cytological staining reagents. For example, a sputum sample is applied to a microscope slide. After appropriate fixation and lysis, hybridization with the probes of the present invention is carried out in situ. In this manner, fungi could be visualized in a specimen by fluorescently labelling Probe 1416 and examining the slide using a fluorescent microscope.

Example 6 Confirmation of Fungemia Following Culture

Following a standard cultivation step utilizing the Bactec, Roche Septi-Chek, or DuPont Isolator, a colony or liquid culture is tested for fungal presence employing Probes 1418 and 1542 in the procedures described in Example 2. Of good advantage is that pure culture is not necessary.

It will be readily appreciated by those skilled in the art that various modifications to the procedures or probes set forth herein may be made without departing from either the spirit or scope of the present invention. In particular, when modifications of the probes such as by deleting one or two end nucleotides with accompanying adjustments in hybridization conditions are to be deemed equivalent.

TABLE 1

DOTBLOT HYBRIDIZATION DATA

NAME/STRAIN		PROBE							
		1417 60degC	1417 65degC	1418 65degC	1415 65degC	1416 65degC	IG707 65degC	1542 60degC	1545 60degC
<i>Alternaria alternata</i>	13963	+	+	+	-	-	-	+	+
<i>Agaricus brunnescens</i>	n5829	+	+	+	+	-	-	+	+
<i>Aspergillus flavus</i>	10124	+	+	+	+	-	-	+	+
<i>Aspergillus fumigatus</i>	36607	+	+	+	+	-	-	+	+
<i>Aspergillus nidulans</i>	10074	+	+	+	-	-	-	+	+
<i>Aspergillus niger</i>	16888	+	+	+	+	-	-	+	+
<i>Aspergillus parasiticus</i>	15517	+	+	+	+	-	-	+	+
<i>Aspergillus terreus</i>	46941	+	+	+	+	-	-	+	+
<i>Aspergillus versicolor</i>	95776	+	+	+	+	-	-	+	+
<i>Blastomyces dermatitidis</i>	60916	+	+	+	+	-	-	+	+
<i>Byssoschlamys fulva</i>	10099	+	+	+	-	-	-	+	+
<i>Candida albicans</i>	11006	+	+	+	-	+	-	+	+
<i>Candida albicans</i>	14053	+	+	+	-	+	-	+	+
<i>Candida albicans</i>	18804	+	+	+	-	+	-	+	+
<i>Candida albicans</i>	24433	+	+	+	-	+	-	+	+
<i>Candida albicans</i>	36232	+	+	+	-	+	-	+	+
<i>Candida albicans</i>	60193	+	+	+	-	+	-	+	+
<i>Candida guilliermondii</i>	6260	+	+	+	-	+	-	+	+
<i>Candida kefyr</i>	4135	+	+	+	-	+	-	+	+
<i>Candida kefyr</i>	46764	+	+	+	-	+	-	+	+
<i>Candida krusei</i>	6258	+	+	+	-	+	-	+	+
<i>Candida lusitanae</i>	42720	+	+	+	-	+	-	+	+
<i>Candida parapsilosis</i>	22019	+	+	+	-	+	-	+	+
<i>Candida rugosa</i>	58964	+	+	+	-	+	-	+	+
<i>Candida tropicalis</i>	750	+	+	+	-	+	-	+	+
<i>Candida tropicalis</i>	13803	+	+	+	-	+	-	+	+
<i>Candida tropicalis</i>	42678	+	+	+	-	+	-	+	+
<i>Candida utilis</i>	9226	+	+	+	-	+	-	+	+
<i>Candida viswanathii</i>	22981	+	+	+	-	+	-	+	+
<i>Chrysosporium keratinophilum</i>	14803	+	+	+	+	-	-	+	+
<i>Cladosporium castellani</i>	24788	+	+	+	+	-	-	+	+

TABLE 1 (CONTINUED)

NAME/STRAIN		PROBE					1545 60degC	1542 60degC	1545 60degC
		1417 60degC	1417 65degC	1418 65degC	1415 65degC	1416 65degC			
Cryptococcus neoformans	14116	+	+	+	+	-	+	+	+
Cryptococcus neoformans	32045	+	+	+	+	-	+	+	+
Cyathus stercoreus	n6473	+	+	+	-	-	+	+	+
Entomophthora virulenta	14207	+	+	+	-	-	+	+	+
Epidermophyton floccosum	52066	+	+	+	+	-	+	+	+
Filobasidiella neoformans	6352	+	+	+	+	-	+	+	+
Fusarium oxysporum	16322	+	+	+	+	-	+	+	+
Hansenula polymorpha	34438	+	+	+	+	+	+	+	+
Histoplasma capsulatum	12700	+	+	+	+	-	+	+	+
Geotrichum candidum	34614	+	+	+	-	-	+	+	+
Lipomyces starkeyi	n11557	+	+	+	-	-	+	+	+
Metschnikowia bicuspidata	22297	+	+	+	-	-	+	+	+
Microsporium racemosum	38556	+	+	+	+	-	+	+	+
Morchella crassipes	18408	+	+	+	+	-	+	+	+
Mucor rouxii	24905	+	+	+	+	-	+	+	+
Neurospora crassa	14692	+	+	+	+	-	+	+	+
Neurospora sitophila	36935	+	+	+	+	-	+	+	+
Paracoccidioides brasiliensis	48093	+	+	+	+	-	+	+	+
Penicillium chrysogenum	10106	+	+	+	+	-	+	+	+
Penicillium digitatum	48113	+	+	+	+	-	+	+	+
Penicillium notatum	9179	+	+	+	+	-	+	+	+
Phycomyces blakesleeanus	n1464	+	+	+	+	+	+	+	+
Pityrosporum ovale	14521	+	+	+	+	-	+	+	+
Pseudallescheria boydii	28169	+	+	+	+	-	+	+	+
Rhizopus oligosporus	22959	+	+	+	+	-	+	+	+
Rhodospiridium toruloides	10788	+	+	+	+	-	+	+	+
Rhodotorula rubra	9449	+	+	+	+	-	+	+	+
Saccharomyces cerevisiae	18824	+	+	+	+	+	+	+	+
Saccharomycodes ludwigii	n12792	+	+	+	+	-	+	+	+
Schizosaccharomyces octosporus	4206	+	+	+	+	-	+	+	+
Sporothrix schenckii	14284	+	+	+	+	-	+	+	+
Taphrina deformans	nT857	+	+	+	+	-	+	+	+

TABLE 1 (CONTINUED)

NAME/STRAIN	1417 60degC	1417 65degC	PROBE 1418 65degC	1415 65degC	1416 65degC	IG707 65degC	1542 60degC	1545 60degC
<i>Torulopsis glabrata</i>	+	+	+	-	+	-	+	+
<i>Tremella mesenterica</i>	+	+	+	-	-	-	+	+
<i>Trichophyton mentagrophytes</i>	+	+	+	-	-	-	+	+
<i>Trichophyton rubrum</i>	+	+	+	-	-	-	+	+
<i>Trichosporon beigeli</i>	+	+	+	+	-	-	+	+
<i>Trichosporon capitatum</i>	+	+	+	-	-	-	+	+
<i>Ustilago maydis</i>	+	+	+	+	-	-	+	+
<i>Verticillium dahliae</i>	+	+	+	+	-	-	+	+
<i>Yarrowia lipolytica</i>	+	+	+	-	-	+	+	+
EXCLUSIVITY								
HUMAN/CaSK1	-	-	-	-	-	-	+	+
<i>Staphylococcus aureus</i>	-	-	-	-	-	-	-	-
<i>Escherichia coli</i>	-	-	-	-	-	-	-	-
Stool RNA	-	nt	-	-	-	-	-	+
Wheat germ RNA	+	nt	+	-	-	-	+	+

TABLE 1 (CONTINUED)

NAME/STRAIN	PROBE					
	1417 60degC	1417 65degC	1418 65degC	1415 65degC	1416 65degC	1545 60degC
<i>Candida albicans</i>						
151-87	+	+	+	-	+	+
184-87	+	+	+	-	+	+
192-87	+	+	+	-	+	+
738-88	+	+	+	-	+	+
784-88	+	+	+	-	+	+
819-88	+	+	+	-	+	+
854-88	+	+	+	-	+	+
864-88	+	+	+	-	+	+
875-88	+	+	+	-	+	+
876-88	+	+	+	-	+	+
889-88	+	+	+	-	+	+
892-88	+	+	+	-	+	+
896-88	+	+	+	-	+	+
901-88	+	+	+	-	+	+
903-88	+	+	+	-	+	+
904-88	+	+	+	-	+	+
917-88	+	+	+	-	+	+
921-88	+	+	+	-	+	+
925-88	+	+	+	-	+	+
926-88	+	+	+	-	+	+
939-88	+	+	+	-	+	+
943-88	+	+	+	-	+	+
946-88	+	+	+	-	+	+
966-88	+	+	+	-	+	+
993-88	+	+	+	-	+	+
161-87	+	+	+	-	+	+
162-87	+	+	+	-	+	+
190-87	+	+	+	-	+	+
203-87	+	+	+	-	+	+
207-87	+	+	+	-	+	+
223-87	+	+	+	-	+	+

TABLE 1 (CONTINUED)

NAME/STRAIN	PROBE					1542 60degC	1545 60degC
	1417 60degC	1417 65degC	1418 65degC	1415 65degC	1416 65degC	1417 65degC	1418 65degC
<i>Candida albicans</i> (continued)							
227-87	+	+	+	-	+	-	+
258-87	+	+	+	-	+	-	+
262-87	+	+	+	-	+	-	+
266-87	+	+	+	-	+	-	+
291-87	+	+	+	-	+	-	+
296-87	+	+	+	-	+	-	+
307-87	+	+	+	-	+	-	+
308-87	+	+	+	-	+	-	+
326-87	+	+	+	-	+	-	+
342-87	+	+	+	-	+	-	+
662-87	+	+	+	-	+	-	+
996-87	+	+	+	-	+	-	+
984-88	+	+	+	-	+	-	+
100888	+	+	+	-	+	-	+
101888	+	+	+	-	+	-	+
<i>Candida guilliermondii</i>							
105586	+	+	+	-	+	-	+
350-87	+	+	+	-	+	-	+
715-88	+	+	+	-	+	-	+
974-88	+	+	+	-	+	-	+
<i>Candida krusei</i>							
46-87	+	+	+	-	+	-	+
528-87	+	+	+	-	+	-	+
842-88	+	+	+	-	+	-	+
939-88	+	+	+	-	+	-	+
<i>Candida (Yarrowia) lipolytica</i>							
056584	+	+	+	-	+	-	+
103486	+	+	+	+	+	+	+
125085	+	+	+	-	+	-	+
453-87	+	+	+	-	+	-	+

TABLE 1 (CONTINUED)

NAME/STRAIN

NAME/STRAIN	PROBE					
	1417 60degC	1417 65degC	1418 65degC	1415 65degC	1416 65degC	1542 60degC
<i>Candida lusitanae</i>						
121585	+	+	+	-	+	+
121585	+	+	+	-	+	+
403-87	+	+	+	-	+	+
964-88	+	+	+	-	+	+
<i>Candida parapsilosis</i>						
175-87	+	+	+	-	+	+
176-87	+	+	+	-	+	+
491-87	+	+	+	-	+	+
492-87	+	+	+	-	+	+
746-88	+	+	+	-	+	+
754-88	+	+	+	-	+	+
828-88	+	+	+	-	+	+
951-88	+	+	+	-	+	+
<i>Candida (kefyr) pseudotropicalis</i>						
091486	+	+	+	-	+	+
100188	+	+	+	-	+	+
102886	+	+	+	-	+	+
999-88	+	+	+	-	+	+
<i>Candida tropicalis</i>						
484-87	+	+	+	-	+	+
784-88	+	+	+	-	+	+
802-88	+	+	+	-	+	+
846-88	+	+	+	-	+	+
997-88	+	+	+	-	+	+
999-88	+	+	+	-	+	+
150-87	+	+	+	-	+	+
210-87	+	+	+	-	+	+
224-87	+	+	+	-	+	+
319-87	+	+	+	-	+	+
573-87	+	+	+	-	+	+

TABLE 1 (CONTINUED)

NAME/STRAIN	PROBE					
	1417 60degC	1417 65degC	1418 65degC	1415 65degC	1416 65degC	1542 60degC
<i>Torulopsis glabrata</i>						
233-87	+	+	+	-	+	+
260-87	+	+	+	-	+	+
275-87	+	+	+	-	+	+
288-87	+	+	+	-	+	+
334-87	+	+	+	-	+	+
359-87	+	+	+	-	+	+
373-87	+	+	+	-	+	+
506-87	+	+	+	-	+	+
562-87	+	+	+	-	+	+
573-87	+	+	+	-	+	+
701-87	+	+	+	-	+	+
901-88	+	+	+	-	+	+
903-88	+	+	+	-	+	+
<i>Cryptococcus albidus</i>						
83-0085	+	nt	+	nt	nt	nt
85-0707	+	nt	+	nt	nt	nt
85-0808	+	nt	+	nt	nt	nt
85-1452	+	nt	+	nt	nt	nt
88-1047	+	nt	+	nt	nt	nt
<i>Cryptococcus laurentii</i>						
82-0600	+	nt	+	nt	nt	nt
87-0657	+	nt	+	nt	nt	nt
88-0010	+	nt	+	nt	nt	nt
<i>Cryptococcus neoformans A</i>						
151	+	nt	+	nt	nt	nt
159	+	nt	+	nt	nt	nt
160	+	nt	+	nt	nt	nt
161	+	nt	+	nt	nt	nt

TABLE 1 (CONTINUED)

NAME/STRAIN	1417 60degC	1417 65degC	PROBE 1418 65degC	1415 65degC	1416 65degC	IG707 65degC	1542 60degC	1545 60degC
Cryptococcus neoformans B								
182	+	nt	+	nt	nt	nt	nt	nt
184	+	nt	+	nt	nt	nt	nt	nt
B3174a	+	nt	+	nt	nt	nt	nt	nt
B3268b	+	nt	+	nt	nt	nt	nt	nt
B3271a	+	nt	+	nt	nt	nt	nt	nt
Cryptococcus neoformans C								
298	+	nt	+	nt	nt	nt	nt	nt
B3185a	+	nt	+	nt	nt	nt	nt	nt
B3186a	+	nt	+	nt	nt	nt	nt	nt
B3267b	+	nt	+	nt	nt	nt	nt	nt
CP110	+	nt	+	nt	nt	nt	nt	nt
Cryptococcus neoformans D								
161	+	nt	+	nt	nt	nt	nt	nt
165C	+	nt	+	nt	nt	nt	nt	nt
166	+	nt	+	nt	nt	nt	nt	nt
167	+	nt	+	nt	nt	nt	nt	nt
168	+	nt	+	nt	nt	nt	nt	nt

DOTBLOT HYBRIDIZATION DATA

NAME/STRAIN	PROBE							
	1859 50degC	1859 60degC	1860 50degC	1860 60degC	1858 50degC	1858 60degC	1857 50degC	1857 60degC
<i>Alternaria alternata</i>	13963	-	+	+	-	-	+	+
<i>Agaricus brunnescens</i>	n5829	-	+	+	-	-	+	+
<i>Aspergillus flavus</i>	10124	-	+	+	-	-	+	+
<i>Aspergillus fumigatus</i>	36607	-	+	+	-	-	+	+
<i>Aspergillus nidulans</i>	10074	-	+	+	-	-	+	+
<i>Aspergillus niger</i>	16888	-	+	+	-	-	+	+
<i>Aspergillus parasiticus</i>	15517	-	+	+	-	-	+	+
<i>Aspergillus terreus</i>	46941	-	+	+	-	-	+	+
<i>Aspergillus versicolor</i>	95776	-	+	+	-	-	+	+
<i>Blastomyces dermatitidis</i>	60916	-	+	+	-	-	+	+
<i>Byssoschlamys fulva</i>	10099	-	+	+	-	-	+	+
<i>Candida albicans</i>	11006	-	+	+	-	-	+	+
<i>Candida albicans</i>	14053	-	+	+	-	-	+	+
<i>Candida albicans</i>	18804	-	+	+	-	-	+	+
<i>Candida albicans</i>	24433	-	+	+	-	-	+	+
<i>Candida albicans</i>	36232	-	+	+	-	-	+	+
<i>Candida albicans</i>	60193	-	+	+	-	-	+	+
<i>Candida guilliermondii</i>	6260	-	+	+	-	-	+	+
<i>Candida kefyr</i>	4135	-	+	+	-	-	+	+
<i>Candida kefyr</i>	46764	-	+	+	-	-	+	+
<i>Candida krusei</i>	6258	-	+	+	-	-	+	+
<i>Candida lusitanae</i>	42720	-	+	+	-	-	+	+
<i>Candida parapsilosis</i>	22019	-	+	+	-	-	+	+
<i>Candida rugosa</i>	58964	-	+	+	-	-	+	+
<i>Candida tropicalis</i>	750	-	+	+	-	-	+	+
<i>Candida tropicalis</i>	13803	-	+	+	-	-	+	+
<i>Candida tropicalis</i>	42678	-	+	+	-	-	+	+
<i>Candida utilis</i>	9226	-	+	+	-	-	+	+
<i>Candida viswanathii</i>	22981	-	+	+	-	-	+	+
<i>Chrysosporium keratinophilum</i>	14803	-	+	+	-	-	+	+
<i>Cladosporium castellani</i>	24788	+	+	+	+	+	+	+

TABLE 2 (CONTINUED)

NAME/STRAIN	PROBE					
	1859 50degC	1859 60degC	1860 50degC	1860 60degC	1858 50degC	1858 60degC
Cryptococcus neoformans 14116	-	-	+	+	-	-
Cryptococcus neoformans 32045	-	-	+	+	+-	-
Cyathus stercoreus n6473	-	-	+	+	+	+-
Entomophthora virulenta 14207	-	-	+	+	+	+
Epidermophyton floccosum 52066	-	-	+	+	+	+
Filobasidiella neoformans 6352	-	-	+	+	+-	+
Fusarium oxysporum 16322	-	-	+	+	+	-
Geotrichum candidum 34614	-	-	+	+	+	+
Hansenula polymorpha 34438	-	-	+	+	+	+
Histoplasma capsulatum 12700	-	-	+	+	+	+
Lipomyces starkeyi n11557	-	-	+	+	+	+
Metschnikowia bicuspidata 22297	-	-	+	+	+	+
Microsporium racemosum 38556	-	-	+	+	+	+
Mucor rouxii 24905	+-	-	+	+	+-	+
Neurospora crassa 14692	-	-	+	+	+	+
Neurospora sitophila 36935	-	-	+	+	+	+
Paracoccidioides brasiliensis 48093	-	-	+	+	+	+
Penicillium chrysogenum 10106	-	-	+	+	+	+
Penicillium digitatum 48113	-	-	+	+	+	+
Penicillium notatum 9179	-	-	+	+	+	+
Phycomyces blakesleeana n1464	-	-	+	+	+-	+-
Pityrosporum ovale 14521	-	-	+	+	-	-
Pseudallescheria boydii 28169	-	-	+	+	-	-
Rhizopus oligosporus 22959	+-	-	+	+	+	+-
Rhodospiridium toruloides 10788	-	-	+	+	+-	-
Rhodotorula rubra 9449	-	-	+	+	-	-
Saccharomyces cerevisiae 18824	-	-	+	+	-	-
Saccharomycodes ludwigii n12792	-	-	+	+	+	+
Schizosaccharomyces octosporus 4206	-	-	+	+	+	+
14284	-	-	+	+	+	+
Sporothrix schenckii nT857	-	-	+	+	+	+
Taphrina deformans 2001	-	-	+	+	+-	+
Torulopsis glabrata	-	-	+	+	-	-

TABLE 2 (CONTINUED)

NAME/STRAIN	PROBE							
	1859 50degC	1859 60degC	1860 50degC	1860 60degC	1858 50degC	1858 60degC	1857 50degC	1857 60degC
<i>Tremella mesenterica</i> 42219	-	-	+	+-	-	-	+	+-
<i>Trichophyton mentagrophytes</i> 28185	-	-	+	+	-	-	+	+
<i>Trichophyton rubrum</i> 28188	-	-	+	+	-	-	+	+
<i>Trichosporon beigeli</i> 28592	-	-	+	+-	-	-	+	+-
<i>Trichosporon capitatum</i> 10663	-	-	+-	-	-	-	+	+
<i>Ustilago maydis</i> 11402	-	-	+	+-	-	-	+	+-
<i>Verticillium dahliae</i> 16535	-	-	+	+	-	-	+	+
<i>Yarrowia lipolytica</i> 18942	+	+-	+	+	-	-	+	+
EXCLUSIVITY								
HUMAN/CaSKI	-	-	+-	-	-	-	+-	-
<i>Staphylococcus aureus</i> GT204	-	-	-	-	-	-	-	-
<i>Escherichia coli</i> 12036	-	-	-	-	-	-	-	-
Stool RNA	-	-	-	-	-	-	-	-
Wheat germ RNA	-	-	+-	+-	-	-	+-	-

TABLE 2 (CONTINUED)

NAME/STRAIN	PROBE					
	1859 50degC	1859 60degC	1860 50degC	1860 60degC	1858 50degC	1857 60degC
<i>Candida albicans</i>						
190-87	-	-	+	+	-	+
266-87	-	-	+	+	-	+
<i>Candida guilliermondii</i>						
350-87	-	-	+	+	-	+
715-88	-	-	+	+	-	+
<i>Candida krusei</i>						
842-88	-	-	+	+	-	+
939-88	-	-	+	+	-	+
<i>Candida (Yarrowia) lipolytica</i>						
056584	+	+	+	+	-	+
125085	+	+	+	+	-	+
<i>Candida lusitanae</i>						
121585	-	-	+	+	-	+
121685	-	-	+	+	-	+
<i>Candida parapsilosis</i>						
176-87	-	-	+	+	-	+
754-88	-	-	+	+	-	+
<i>Candida (kefir) pseudotropicalis</i>						
100188	-	-	+	+	+	+
999-88	-	-	+	+	+	+
<i>Candida tropicalis</i>						
150-87	-	-	+	+	-	+
210-87	-	-	+	+	-	+

TABLE 2 (CONTINUED)

NAME/STRAIN	PROBE					
	1859 50degC	1859 60degC	1860 50degC	1860 60degC	1858 50degC	1857 60degC
<i>Torulopsis glabrata</i>						
288-87	-	-	+	+	+	+
334-87	-	-	+	+	+	+
<i>Cryptococcus albidus</i>						
83-0085	-	-	+	+	+	-
85-0707	-	-	+	+	+	-
<i>Cryptococcus laurentii</i>						
82-0600	-	-	+	+	+	+
87-0657	-	-	+	+	+	+
<i>Cryptococcus neoformans</i> A						
151	-	-	+	+	+	-
162	-	-	+	+	+	-
<i>Cryptococcus neoformans</i> B						
184	-	-	+	+	+	-
B3174a	-	-	+	+	+	-
<i>Cryptococcus neoformans</i> C						
298	-	-	+	+	+	-
B3185a	-	-	+	+	+	-
<i>Cryptococcus neoformans</i> D						
161	-	-	+	+	+	-
166	-	-	+	+	+	-

TABLE 2 (CONTINUED)

NAME/STRAIN	1859 50degC	1859 60degC	PROBE 1860 50degC	1860 60degC	1858 50degC	1858 60degC	1857 50degC	1857 60degC
ADDITIONAL ASPERGILLI AND PENICILLIUM								
Aspergillus clavatus 9192	-	-	+	+	-	-	+	+
Aspergillus niger 10535	-	-	+	+	-	-	+	+
Aspergillus niger 16404	-	-	+	+	-	-	+	+
Aspergillus niger 16880	-	-	+	+	-	-	+	+
Aspergillus oryzae 14895	-	-	+	+	-	-	+	+
Aspergillus repens 48521	-	-	+	+	-	-	+	+
Penicillium chrysogenum 10002	-	-	+	+	-	-	+	+
Penicillium expansum 7861	-	-	+	+	-	-	+	+
Penicillium italicum 48114	-	-	+	+	-	-	+	+
Penicillium roquefortii 10110	-	-	+	+	-	-	+	+
ADDITIONAL BLASTOMYCES AND HISTOPLASMAS								
Blastomyces dermatitidis 064571	-	-	+	+	-	-	+	+
Blastomyces dermatitidis 82-0741	-	-	+	+	-	-	+	+
Blastomyces dermatitidis CDC B44	-	-	+	+	-	-	+	+
Blastomyces dermatitidis NYS48	-	-	+	+	-	-	+	+
Histoplasma capsulatum CDC A28	-	-	+	+	-	-	+	+
Histoplasma capsulatum NYS 211	-	-	+	+	-	-	+	+
Histoplasma capsulatum NYS 214	-	-	+	+	-	-	+	+
Histoplasma capsulatum NYS 215	-	-	+	+	-	-	+	+
Histoplasma capsulatum NYS 216	-	-	+	+	-	-	+	+
Histoplasma capsulatum NYS 232	-	-	+	+	-	-	+	+

TABLE 3

DOTBLOT HYBRIDIZATION DATA

NAME/STRAIN		PROBE							
		1813 50degC	1813 60degC	1814 50degC	1814 60degC	1812 50degC	1812 60degC	1816 50degC	1816 60degC
<i>Alternaria alternata</i>	13963	+	-	+-	-	+	+-	+-	-
<i>Agaricus brunnescens</i>	n5829	+	+-	+	+	+	+	+	-
<i>Aspergillus flavus</i>	10124	+	-	+	+	+	+-	+-	-
<i>Aspergillus fumigatus</i>	36607	+	-	+	+	+	+-	+-	-
<i>Aspergillus nidulans</i>	10074	+	-	+	+	+	+-	+-	-
<i>Aspergillus niger</i>	16888	+	-	+	+-	+	+-	+-	-
<i>Aspergillus parasiticus</i>	15517	+	-	+	+	+	+-	+-	-
<i>Aspergillus terreus</i>	46941	+	-	+	+	+	+-	+-	-
<i>Aspergillus versicolor</i>	95776	+	-	+	+-	+	+-	+-	-
<i>Blastomyces dermatitidis</i>	60916	+	-	+	+	+	+-	+	-
<i>Byssoschlamys fulva</i>	10099	+	-	+	+	+	+-	+-	-
<i>Candida albicans</i>	11006	+	-	+	+	+	+-	+	-
<i>Candida albicans</i>	14053	+	-	+	+	+	+-	+	-
<i>Candida albicans</i>	18804	+	-	+	+	+	+-	+	-
<i>Candida albicans</i>	24433	+	-	+	+	+	+-	+	-
<i>Candida albicans</i>	36232	+	-	+	+	+	+-	+	-
<i>Candida albicans</i>	60193	+	-	+	+	+	+-	+	-
<i>Candida guilliermondii</i>	6260	+	-	+	+	+	+-	+	-
<i>Candida kefir</i>	4135	+	-	+	+	+	+-	+	-
<i>Candida kefir</i>	46764	+	-	+	+	+	+-	+	-
<i>Candida krusei</i>	6258	+	-	+	-	+	+-	+	-
<i>Candida lusitanae</i>	42720	+	-	+	-	+	+-	+	-
<i>Candida parapsilosis</i>	22019	+	-	+	+	+	+-	+	-
<i>Candida rugosa</i>	58964	+	+-	+-	-	+	+-	+-	-
<i>Candida tropicalis</i>	750	+	-	+	+	+	+-	+	-
<i>Candida tropicalis</i>	13803	+	-	+	+	+	+-	+	-
<i>Candida tropicalis</i>	42678	+	-	+	+	+	+-	+	-
<i>Candida utilis</i>	9226	+	-	+	+-	+	+-	+	-
<i>Candida viswanathii</i>	22981	+	-	+	+	+	+-	+	-
<i>Chrysosporium keratinophilum</i>	14803	+	-	+	+	+	+-	+	-
<i>Cladosporium castellani</i>	24788	+	-	+	+	+	+-	+	-

TABLE 3 (CONTINUED)

NAME/STRAIN	PROBE							
	1813 50degC	1813 60degC	1814 50degC	1814 60degC	1812 50degC	1812 60degC	1816 50degC	1816 60degC
Cryptococcus neoformans	+	-	+	+	+	++	-	-
Cryptococcus neoformans	+	-	+	+	+	++	-	-
Cyathus stercoreus	+	-	+	+	+	+	+	-
Entomophthora virulenta	+	-	+	+	+	++	+	-
Epidermophyton floccosum	+	-	+	+	+	++	+	-
Filobasidiella neoformans	+	-	+	+	+	++	+	+
Fusarium oxysporum	+	-	+	+	+	++	+	-
Geotrichum candidum	+	-	+	+	+	++	+	-
Hansenula polymorpha	+	-	+	+	+	++	+	-
Histoplasma capsulatum	+	-	+	+	+	++	+	-
Lipomyces starkeyi	+	-	+	+	+	++	+	-
Metschnikowia bicuspidata	+	-	+	+	+	+	+	-
Microsporium racemosum	+	-	+	+	+	++	+	-
Mucor rouxii	+	-	+	+	+	-	+	-
Neurospora crassa	+	-	+	+	+	++	+	+
Neurospora sitophila	+	-	+	+	+	++	+	+
Paracoccidioides brasiliensis	+	-	+	+	+	++	+	-
Penicillium chrysogenum	+	-	+	+	+	++	+	-
Penicillium digitatum	+	-	+	+	+	-	+	-
Penicillium notatum	+	-	+	+	+	++	+	-
Phycomyces blakesleeana	+	-	+	+	+	++	+	-
Pityrosporum ovale	+	-	+	+	+	+	+	+
Pseudallescheria boydii	+	-	+	+	+	++	+	+
Rhizopus oligosporus	+	-	+	+	+	-	+	-
Rhodospiridium toruloides	+	-	+	+	+	++	+	-
Rhodotorula rubra	+	-	+	+	+	+	+	-
Saccharomyces cerevisiae	+	-	+	+	+	++	+	-
Saccharomyces ludwigii	+	-	+	+	+	++	+	-
Schizosaccharomyces octosporus	+	-	+	+	+	++	+	-
Sporothrix schenckii	+	-	+	+	+	++	+	-
Taphrina deformans	+	-	+	+	+	++	+	-
Torulopsis glabrata	+	-	+	+	+	++	+	-

TABLE 3 (CONTINUED)

NAME/STRAIN	PROBE							
	1813 50degC	1813 60degC	1814 50degC	1814 60degC	1812 50degC	1812 60degC	1816 50degC	1816 60degC
<i>Tremella mesenterica</i>	-	-	+	+	+	-	+	-
<i>Trichophyton mentagrophytes</i>	+	-	+	+	+	+	+	-
<i>Trichophyton rubrum</i>	+	-	+	+	+	+	+	-
<i>Trichosporon beigellii</i>	+	-	+	+	+	+	+	-
<i>Trichosporon capitatum</i>	-	-	+	+	+	-	+	-
<i>Ustilago maydis</i>	-	-	+	+	+	+	+	-
<i>Verticillium dahliae</i>	+	-	+	+	+	+	+	+
<i>Yarrowia lipolytica</i>	+	-	+	+	+	-	+	-
42219								
28185								
28188								
28592								
10663								
11402								
16535								
18942								
EXCLUSIVITY								
HUMAN/CASKI								
<i>Staphylococcus aureus</i>								
<i>Escherichia coli</i>								
Stool RNA								
Wheat germ RNA								
GT204								
12036								

TABLE 3 (CONTINUED)

NAME/STRAIN	PROBE					
	1813 50degC	1813 60degC	1814 50degC	1814 60degC	1812 50degC	1812 60degC
<i>Candida albicans</i>						
190-87	+	-	+	+	+	+
266-87	+	-	+	+	+	+
<i>Candida guilliermondii</i>						
350-87	+	-	+	+	+	+
715-88	+	-	+	+	+	+
<i>Candida krusei</i>						
842-88	+	-	+	+	+	+
939-88	+	-	+	+	+	+
<i>Candida (Yarrowia) lipolytica</i>						
056584	+-	-	-	-	-	-
125085	+-	-	-	-	-	-
<i>Candida lusitanae</i>						
121585	+	-	+	+	+	+
121685	+	-	+	+	+	+
<i>Candida parapsilosis</i>						
176-87	+	-	+	+	+	+
754-88	+	-	+	+	+	+
<i>Candida (kefyr) pseudotropicalis</i>						
100188	+	-	+	+	+	+
999-88	+	-	+	+	+	+
<i>Candida tropicalis</i>						
150-87	+	-	+	+	+	+
210-87	+	-	+	+	+	+

TABLE 3 (CONTINUED)

NAME/STRAIN	PROBE					
	1813 50degC	1813 60degC	1814 50degC	1814 60degC	1812 50degC	1816 60degC
<i>Torulopsis glabrata</i>						
288-87	+	-	+	+	+-	-
334-87	+	-	+	+	+-	-
<i>Cryptococcus albidus</i>						
83-0085	+	-	+	+	+-	-
85-0707	+	-	+	+	+-	-
<i>Cryptococcus laurentii</i>						
82-0600	+-	-	+	+	+-	-
87-0657	+-	-	+	+	+-	-
<i>Cryptococcus neoformans</i> A						
151	+	-	+	+	+-	-
162	+	-	+	+	+-	-
<i>Cryptococcus neoformans</i> B						
184	+	-	+	+	+-	-
B3174a	+	-	+	+	+-	-
<i>C-ryptococcus neoformans</i> C						
298	+	-	+	+	+-	-
B3185a	+	-	+	+	+-	-
<i>Cryptococcus neoformans</i> D						
161	+	-	+	+	+-	-
166	+	-	+	+	+-	-

TABLE 3 (CONTINUED)

NAME/STRAIN	PROBE							
	1813 50degC	1813 60degC	1814 50degC	1814 60degC	1812 50degC	1812 60degC	1816 50degC	1816 60degC
ADDITIONAL ASPERGILLI AND PENICILLIUM								
Aspergillus clavatus	+	-	+	+	+	+-	nt	nt
Aspergillus niger	+	-	+	+	+	+-	nt	nt
Aspergillus niger	+	-	+	+	+	+-	nt	nt
Aspergillus niger	+	-	+	+	+	+-	nt	nt
Aspergillus oryzae	+	-	+	+	+	+-	nt	nt
Aspergillus repens	+	-	+	+	+	+-	nt	nt
Penicillium chrysogenum	+	-	+	+	+	+-	nt	nt
Penicillium expansum	+	-	+	+	+	+-	nt	nt
Penicillium italicum	+	-	+	+	+	+-	nt	nt
Penicillium roquefortii	+	-	+	+	+	+-	nt	nt
ADDITIONAL BLASTOMYCES AND HISTOPLASMAS								
Blastomyces dermatitidis 064571	+	-	+	+	+	+-	nt	nt
Blastomyces dermatitidis 82-0741	+	-	+	+	+	+-	nt	nt
Blastomyces dermatitidis CDC B44	+	-	+	+	+	+-	nt	nt
Blastomyces dermatitidis NYS48	+	-	+	+	+	+-	nt	nt
Histoplasma capsulatum CDC A28	+	-	+	+	+	+-	nt	nt
Histoplasma capsulatum NYS 211	+	-	+	+	+	+-	nt	nt
Histoplasma capsulatum NYS 214	+	-	+	+	+	+-	nt	nt
Histoplasma capsulatum NYS 215	+	-	+	+	+	+-	nt	nt
Histoplasma capsulatum NYS 216	+	-	+	+	+	+-	nt	nt
Histoplasma capsulatum NYS 232	+	-	+	+	+	+-	nt	nt

SEQUENCE LISTING

SEQ ID NO.1

SEQUENCE TYPE: Nucleotide
 SEQUENCE LENGTH: 27 bases
 STRANDEDNESS: single
 TOPOLOGY: Linear
 MOLECULAR TYPE: Synthetic or biological DNA or RNA
 HYPOTHETICAL?: No
 SENSE OR ANTI-SENSE?: Sense
 FEATURES: Referred to in the specification as probe 1417
 PROPERTIES: Probe for selective detection of fungus.
 SEQUENCE:
 TGTCTGGACC TGGTGAGTTT CCCCCTG

SEQUENCE LISTING

SEQ ID NO.2

SEQUENCE TYPE: Nucleotide
 SEQUENCE LENGTH: 39 bases
 STRANDEDNESS: single
 TOPOLOGY: Linear
 MOLECULAR TYPE: Synthetic or biological DNA or RNA
 HYPOTHETICAL?: No
 SENSE OR ANTI-SENSE?: Sense
 FEATURES: Referred to in the specification as probe 1418
 PROPERTIES: Probe for selective detection of fungus.
 SEQUENCE:
 TGTCTGGACC TGGTGAGTTT CCCCCTGTTG
 AGTCAAATT

SEQUENCE LISTING

SEQ ID NO.3

SEQUENCE TYPE: Nucleotide
 SEQUENCE LENGTH: 37 bases
 STRANDEDNESS: single
 TOPOLOGY: Linear
 MOLECULAR TYPE: Synthetic or biological DNA or RNA
 HYPOTHETICAL?: No
 SENSE OR ANTI-SENSE?: Sense
 FEATURES: Referred to in the specification as probe 1415
 PROPERTIES: Probe for selective detection of fungus.
 SEQUENCE:

TCCTCGTTAA GGGATTTAAA TTGTACTCAT
 TCCAATT

5 SEQUENCE LISTING

SEQ ID NO.4

SEQUENCE TYPE: Nucleotide
 SEQUENCE LENGTH: 37 bases
 STRANDEDNESS: single
 TOPOLOGY: Linear
 MOLECULAR TYPE: Synthetic or biological DNA or RNA
 HYPOTHETICAL?: No
 SENSE OR ANTI-SENSE?: Sense
 FEATURES: Referred to in the specification as probe 1416
 PROPERTIES: Probe for selective detection of fungus.
 SEQUENCE:
 TCCTCGTTAA GGTATTTACA TTGTACTCAT
 TCCAATT

25 SEQUENCE LISTING

SEQ ID NO.5

SEQUENCE TYPE: Nucleotide
 SEQUENCE LENGTH: 37 bases
 STRANDEDNESS: single
 TOPOLOGY: Linear
 MOLECULAR TYPE: Synthetic or biological DNA or RNA
 HYPOTHETICAL?: No
 SENSE OR ANTI-SENSE?: Sense
 FEATURES: Referred to in the specification as probe IG707
 PROPERTIES: Probe for selective detection of fungus.
 SEQUENCE:
 TCCTCGTTAA GGTGTTTAAA TTGTACTCAT
 TCCAATT

45 SEQUENCE LISTING

SEQ ID NO.6

SEQUENCE TYPE: Nucleotide
 SEQUENCE LENGTH: 26 bases
 STRANDEDNESS: single
 TOPOLOGY: Linear
 MOLECULAR TYPE: Synthetic or biological DNA or RNA
 HYPOTHETICAL?: No
 SENSE OR ANTI-SENSE?: Sense
 FEATURES: Referred to in the specification as probe 1542

PROPERTIES: Probe for selective detection of fungus.

SEQUENCE:

AACTAAGAAC GGCCATGCAC CACCAT

SEQUENCE LISTING

SEQ ID NO.7

SEQUENCE TYPE: Nucleotide

SEQUENCE LENGTH: 40 bases

STRANDEDNESS: single

TOPOLOGY: Linear

MOLECULAR TYPE: Synthetic or biological DNA or RNA

HYPOTHETICAL?: No

SENSE OR ANTI-SENSE?: Sense

FEATURES: Referred to in the specification as probe 1545

PROPERTIES: Probe for selective detection of fungus.

SEQUENCE:

TGGTGCCCTT CCGTCAATTT CTTTAAGTTT
CAGCCTTGCG

SEQUENCE LISTING

SEQ ID NO.8

SEQUENCE TYPE: Nucleotide

SEQUENCE LENGTH: 48 bases

STRANDEDNESS: single

TOPOLOGY: Linear

MOLECULAR TYPE: Synthetic or biological DNA or RNA

HYPOTHETICAL?: No

SENSE OR ANTI-SENSE?: Sense

FEATURES: Referred to in the specification as probe 1814

PROPERTIES: Probe for selective detection of fungus.

SEQUENCE:

TCGCTGGCGC AAGGCCATGC GATTCGAGAG
GTTATTATGA ATCATCAG

SEQUENCE LISTING

SEQ ID NO.9

SEQUENCE TYPE: Nucleotide

SEQUENCE LENGTH: 30 bases

STRANDEDNESS: single

TOPOLOGY: Linear

MOLECULAR TYPE: Synthetic or biological DNA or RNA

HYPOTHETICAL?: No

SENSE OR ANTI-SENSE?: Sense

FEATURES: Referred to in the specification as

probe 1816

PROPERTIES: Probe for selective detection of fungus.

SEQUENCE:

5 CAAGCTGATG ACTTGTGCTT ACTAGGGATT

SEQUENCE LISTING

SEQ ID NO.10

SEQUENCE TYPE: Nucleotide

SEQUENCE LENGTH: 37 bases

STRANDEDNESS: single

TOPOLOGY: Linear

15 MOLECULAR TYPE: Synthetic or biological DNA or RNA

HYPOTHETICAL?: No

SENSE OR ANTI-SENSE?: Sense

FEATURES: Referred to in the specification as probe 1857

20 PROPERTIES: Probe for selective detection of fungus.

SEQUENCE:

25 TCGGCATAGT TTGTGGTTAA GACTACGACG
GTATCTT

SEQUENCE LISTING

SEQ ID NO.11

SEQUENCE TYPE: Nucleotide

SEQUENCE LENGTH: 25 bases

STRANDEDNESS: single

TOPOLOGY: Linear

35 MOLECULAR TYPE: Synthetic or biological DNA or RNA

HYPOTHETICAL?: No

SENSE OR ANTI-SENSE?: Sense

FEATURES: Referred to in the specification as probe 1813

40 PROPERTIES: Probe for selective detection of fungus.

SEQUENCE:

45 AAATGCTTTC GCAGTAGTTG GTCTT

SEQUENCE LISTING

SEQ ID NO.12

50 SEQUENCE TYPE: Nucleotide

SEQUENCE LENGTH: 49 bases

STRANDEDNESS: single

TOPOLOGY: Linear

55 MOLECULAR TYPE: Synthetic or biological DNA or RNA

HYPOTHETICAL?: No

SENSE OR ANTI-SENSE?: Sense

FEATURES: Referred to in the specification as

probe 1860

PROPERTIES: Probe for selective detection of fungus.

SEQUENCE:

AAATGCTTTC GCAGTAGTTG GTCTTCGGTA 5
AATCCAAGAA TTTCACCTT

SEQUENCE LISTING

SEQ ID NO.13

SEQUENCE TYPE: Nucleotide

SEQUENCE LENGTH: 28 bases

STRANDEDNESS: single

TOPOLOGY: Linear

MOLECULAR TYPE: Synthetic or biological DNA or RNA

HYPOTHETICAL?: No

SENSE OR ANTI-SENSE?: Sense

FEATURES: Referred to in the specification as probe 1812

PROPERTIES: Probe for selective detection of fungus.

SEQUENCE:

ACGTCCTATT TTATTATTCC ATGCTAAT

SEQUENCE LISTING

SEQ ID NO.14

SEQUENCE TYPE: Nucleotide

SEQUENCE LENGTH: 29 bases

STRANDEDNESS: single

TOPOLOGY: Linear

MOLECULAR TYPE: Synthetic or biological DNA or RNA

HYPOTHETICAL?: No

SENSE OR ANTI-SENSE?: Sense

FEATURES: Referred to in the specification as probe 1858

PROPERTIES: Probe for selective detection of fungus.

SEQUENCE:

AAGTCATATT TCATTATTCC ATGCTAACT

SEQUENCE LISTING

SEQ ID NO. 15

SEQUENCE TYPE: Nucleotide

SEQUENCE LENGTH: 28 bases

STRANDEDNESS: single

TOPOLOGY: Linear

MOLECULAR TYPE: Synthetic or biological DNA or RNA

HYPOTHETICAL?: No

SENSE OR ANTI-SENSE?: Sense

FEATURES: Referred to in the specification as

probe 1859

PROPERTIES: Probe for selective detection of fungus.

SEQUENCE:

TCGTCGAGTT ATGTTATTCC ATGCAAAT

SEQUENCE LISTING

SEQ ID NO. 16

SEQUENCE TYPE: Nucleotide

SEQUENCE LENGTH: 35 bases

STRANDEDNESS: single

TOPOLOGY: Linear

MOLECULAR TYPE: Synthetic or biological DNA or RNA

HYPOTHETICAL?: No

SENSE OR ANTI-SENSE?: Sense

FEATURES: Referred to in the specification as probe/primer 936

PROPERTIES: Probe for selective detection of fungus.

SEQUENCE:

CCGAATTCGT CGACAACCTG GTTGATCCTG 25
CCAGT

SEQUENCE LISTING

SEQ ID NO. 17

SEQUENCE TYPE: Nucleotide

SEQUENCE LENGTH: 39 bases

STRANDEDNESS: single

TOPOLOGY: Linear

MOLECULAR TYPE: Synthetic or biological DNA or RNA

HYPOTHETICAL?: No

SENSE OR ANTI-SENSE?: Sense

FEATURES: Referred to in the specification as probe/primer 935

PROPERTIES: Probe for selective detection of fungus.

SEQUENCE:

CCCGGGATCC AAGCTTGATC CTTCTGCAGG 45
TTCACCTAC

Claims

50 1. A nucleic acid fragment capable of hybridizing to rRNA or rDNA of fungi but not to rRNA or rDNA of human, bacteria, or wheat.

2. A nucleic acid fragment as claimed in claim 1 wherein said fragment is complementary to at least 90% of a sequence comprising any ten consecutive nucleotides within probes 1417, 1418, 1415, 1416, 1542, 1545, 1G707, 1859, 1860, 1858, 1857, 1812, 1813, 1814, 1816, 936 or 935.

3. A nucleic acid fragment as claimed in claim 1, wherein said fragment is homologous to at least 90% of a sequence comprising any ten consecutive nucleotides within probes 1417, 1418, 1415, 1416, 1542, 1545, IG707, 1859, 1860, 1858, 1857, 1812, 1813, 1814, 1816, 936 or 935.

4. A set of probes comprising at least two nucleic acid fragments, at least one of which comprises probe 1417, 1418, 1415, 1416, 1542, 1545, IG707, 1859, 1860, 1858, 1857, 1812, 1813, 1814, 1816, 936 or 935 or one of their complementary sequences.

5. A nucleic acid fragment as claimed in claim 1 which comprises probe 1417, 1418, 1415, 1416, 1542, 1545, IG707, 1859, 1860, 1858, 1857, 1812, 1813, 1814, 1816, 936 or 935 or one of their complementary sequences.

6. A method for detecting fungal organisms in a sample comprising:

6. A method for detecting fungal organisms in a sample comprising:

a) contacting said sample with at least one nucleic acid fragment under conditions that allow said fragment to hybridize to rRNA or rDNA of said fungal organism if present in said sample, whereby nucleic acid complexes are formed, and wherein said nucleic acid fragment does not hybridize to rRNA or rDNA of non-fungal organisms; and

b) detecting said nucleic acid complexes as an indication of the presence of said fungal organism.

7. A method as claimed in claim 6 wherein said nucleic acid fragment in said contacting step comprises probe 1417, 1418, 1415, 1416, 1542, 1545, IG707, 1859, 1860, 1858, 1857, 1812, 1813, 1814 or 1816 or one of their complementary sequences.

8. A method as claimed in Claim 6 wherein said nucleic acid fragment in said contacting step comprises probe/primer 936 and said detecting step comprises further contacting said sample with a second nucleic acid fragment comprising probe 935, 1417, 1418, 1415, 1416, 1542, 1545, IG707, 1859, 1860, 1858, 1857, 1812, 1813, 1814 or 1816.

9. A method as claimed in Claim 7 further comprising the step of amplifying 18S rRNA gene sequences of said fungal organism by polymerase chain reaction.

Claims for the following Contracting State: ES

1. A process for the preparation of a nucleic acid fragment capable of hybridizing to rRNA or rDNA of fungi but not to rRNA or rDNA of human, bacteria, or wheat, the process comprising coupling successive nucleotides together and/or ligating oligo- and/or poly-nucleotides.

2. A process as claimed in claim 1, wherein said

fragment is complementary to at least 90% of a sequence comprising any ten consecutive nucleotides within probes 1417, 1418, 1415, 1416, 1542, 1545, IG707, 1859, 1860, 1858, 1857, 1812, 1813, 1814, 1816, 936 or 935.

3. A process as claimed in claim 1 wherein said fragment is homologous to at least 90% of a sequence comprising any ten consecutive nucleotides within probes 1417, 1418, 1415, 1416, 1542, 1545, IG707, 1859, 1860, 1858, 1857, 1812, 1813, 1814, 1816, 936 or 935.

4. A process for the preparation of a set of probes comprising at least two nucleic acid fragments, at least one of which comprises probe 1417, 1418, 1415, 1416, 1542, 1545, IG707, 1859, 1860, 1858, 1857, 1812, 1813, 1814, 1816, 936 or 935 or one of their complementary sequences, the process comprising coupling successive nucleotides together and/or ligating oligo- and/or poly-nucleotides.

5. A process as claimed in claim 1, wherein the nucleic acid fragment comprises probe 1417, 1418, 1415, 1416, 1542, 1545, IG707, 1859, 1860, 1858, 1857, 1812, 1813, 1814, 1816, 936 or 935 or one of their complementary sequences.

6. A method for detecting fungal organisms in a sample comprising:

a) contacting said sample with at least one nucleic acid fragment under conditions that allow said fragment to hybridize to rRNA or rDNA of said fungal organism if present in said sample, whereby nucleic acid complexes are formed, and wherein said nucleic acid fragment does not hybridize to rRNA or rDNA of non-fungal organisms; and


b) detecting said nucleic acid complexes as an indication of the presence of said fungal organism.

7. A method as claimed in claim 6 wherein said nucleic acid fragment in said contacting step comprises probe 1417, 1418, 1415, 1416, 1542, 1545, IG707, 1859, 1860, 1858, 1857, 1812, 1813, 1814 or 1816 or one of their complementary sequences.

8. A method as claimed in claim 6 wherein said nucleic acid fragment in said contacting step comprises probe/primer 936 and said detecting step comprises further contacting said sample with a second nucleic acid fragment comprising probe 935, 1417, 1418, 1415, 1416, 1542, 1545, IG707, 1859, 1860, 1858, 1857, 1812, 1813, 1814 or 1816.

9. A method as claimed in claim 7 further comprising the step of amplifying 18S rRNA gene sequences of said fungal organism by polymerase chain reaction.

FIGURE



(19)



Europäisches Patentamt
European Patent Office
Office européen des brevets




(11) Publication number:

0 422 872 A3

(12)

EUROPEAN PATENT APPLICATION(21) Application number: **90311003.9**(51) Int. Cl.⁵: **C12Q 1/68, C07H 21/04**(22) Date of filing: **08.10.90**(30) Priority: **12.10.89 US 420577**(43) Date of publication of application:
17.04.91 Bulletin 91/16(84) Designated Contracting States:
AT BE CH DE DK ES FR GB GR IT LI LU NL SE(88) Date of deferred publication of the search report:
05.01.94 Bulletin 94/01(71) Applicant: **GENE-TRAK SYSTEMS**
31 New York Avenue
Framingham Massachusetts 01701(US)(72) Inventor: **Weisburg, William G.**
3 Jillson Circle**Milford, Massachusetts 01757(US)**Inventor: **Sogin, Mitchell I.****1371 GlenCoe Street****Denver, Colorado 80220(US)**Inventor: **Barns, Susan M.****3 Maple Street****Hopkinton, Massachusetts 01748(US)**Inventor: **Pelletier, Dale A.****58 Walnut Street,****Apartment 3****Somerville, Massachusetts 02143(US)**(74) Representative: **Sheard, Andrew Gregory et al**
Kilburn & Strode
30, John Street
London WC1N 2DD (GB)(54) **Nucleic acid probes and methods for detecting fungi.**

(57) Nucleic acid probes are described for detecting fungi capable of causing fungal septicemia or capable of causing food spoilage. The preferred probes are complementary to ribonucleic acid sequences found in numerous fungi and absent in animal or plant genomes. As such, these probes can detect the rRNA, rDNA, or polymerase chain reaction amplification products from the majority of fungal species. The detection of etiological agents of human fungemia, the clinical diagnosis of this disease and the direct evaluation of food or beverage fungal content utilizing rRNA or rDNA probes is now possible.

FIGURE



European Patent
Office

EUROPEAN SEARCH REPORT

Application Number
EP 90 31 1003.9

DOCUMENTS CONSIDERED TO BE RELEVANT			
Category	Citation of document with indication, where appropriate, of relevant passages	Relevant to claim	CLASSIFICATION OF THE APPLICATION (Int.Cl.5)
D,X	EP-A-0 272 009 (HOGAN ET AL.) * the whole document * * especially page 71-74, claims 241-265; page 139, claim 69.* ---	1-3,6	C12Q1/68 C07H21/04
X Y	EP-A-0 335 633 (INTEGRATED GENETICS INC.) * the whole document * ---	1,6 1-9	
X	EP-A-0 256 843 (CETUS CORP.) * page 5, line 15 * ---	2,3	
E	EP-A-0 421 725 (GENE-TRAK SYST.) * page 6 * ---	2,3	
P,X	WO-A-90 09180 (UNIV. TEXAS SYST) * page 39; claim 8 * ---	2,3	
Y	NUCLEIC ACIDS RESEARCH vol. 17, no. 14, July 1989, OXFORD, GB pages 7843 - 7853 EDWARDS, U. ET AL. 'isolation and direct complete nucleotide determination of entire genes.' * page 7849 - page 7852 * ---	8,9	TECHNICAL FIELDS SEARCHED (Int.Cl.5) C12Q
Y	EP-A-0 314 294 (GENE-TRAK SYSTEMS) * the whole document * ---	1-9	
E	WO-A-91 02092 (GENE-TRAK SYST) * page 11 * * with respect to sequence 1860* ---	2,3	
X	EP-A-0 308 716 (GRUPPO LEPETIT SPA) * page 11; figure 18 * * with respect to sequence 1858* ---	2,3	
-/--			
The present search report has been drawn up for all claims			
Place of search THE HAGUE		Date of completion of the search 4 November 1993	Examiner OSBORNE, H
CATEGORY OF CITED DOCUMENTS			
X : particularly relevant if taken alone Y : particularly relevant if combined with another document of the same category A : technological background O : non-written disclosure P : intermediate document		T : theory or principle underlying the invention E : earlier patent document, but published on, or after the filing date D : document cited in the application L : document cited for other reasons & : member of the same patent family, corresponding document	



European Patent
Office

CLAIMS INCURRING FEES

The present European patent application comprised at the time of filing more than ten claims.

- ☐ All claims fees have been paid within the prescribed time limit. The present European search report has been drawn up for all claims.
- ☐ Only part of the claims fees have been paid within the prescribed time limit. The present European search report has been drawn up for the first ten claims and for those claims for which claims fees have been paid,
namely claims:
- ☐ No claims fees have been paid within the prescribed time limit. The present European search report has been drawn up for the first ten claims.

LACK OF UNITY OF INVENTION

The Search Division considers that the present European patent application does not comply with the requirement of unity of invention and relates to several inventions or groups of inventions,
namely:

See sheet -B-

- ☒ All further search fees have been paid within the fixed time limit. The present European search report has been drawn up for all claims.
- ☐ Only part of the further search fees have been paid within the fixed time limit. The present European search report has been drawn up for those parts of the European patent application which relate to the inventions in respect of which search fees have been paid,
namely claims:
- ☐ None of the further search fees has been paid within the fixed time limit. The present European search report has been drawn up for those parts of the European patent application which relate to the invention first mentioned in the claims,
namely claims:



European Patent
Office

EUROPEAN SEARCH REPORT

Application Number
EP 90 31 1003

DOCUMENTS CONSIDERED TO BE RELEVANT			
Category	Citation of document with indication, where appropriate, of relevant passages	Relevant to claim	CLASSIFICATION OF THE APPLICATION (Int.Cl.5)
X	PROCEEDINGS OF THE NATIONAL ACADEMY OF SCIENCES OF USA. vol. 81, January 1984, WASHINGTON US pages 206 - 210 ANDERSON, R. ET AL 'a pattern of partially homologous recombination in mouse L cells' * figures 2,J6 * * with respect to sequence 1415* ----	2,3	
X	WO-A-87 01132 (KIRIN-AMEGEN) *with respect to sequence 1857* * page 18 * ----	2,3	
X	WO-A-86 00926 (CELLTECH LTD.) * page 20; claim 7 * * with respect to sequence 1859* ----	2,3	
A	WO-A-89 05359 (INTEGRATED GENETICS INC.) ----		
A	WO-A-89 06704 (MICROPROBE CORP.) -----		
The present search report has been drawn up for all claims			TECHNICAL FIELDS SEARCHED (Int.Cl.5)
Place of search THE HAGUE		Date of completion of the search 4 November 1993	Examiner OSBORNE, H
CATEGORY OF CITED DOCUMENTS			
X : particularly relevant if taken alone Y : particularly relevant if combined with another document of the same category A : technological background O : non-written disclosure P : intermediate document		T : theory or principle underlying the invention E : earlier patent document, but published on, or after the filing date D : document cited in the application L : document cited for other reasons & : member of the same patent family, corresponding document	



European Patent
Office

EP 90 31 1003 -B-

LACK OF UNITY OF INVENTION

The Search Division considers that the present European patent application does not comply with the requirement of unity of invention and relates to several inventions or groups of inventions, namely:

1. Claims 1,6 and partially claims 2-5,7-9 where they relate to probes 1417, 1418, 1542, 1545, 936 and 935; probes which are 100 % inclusive for 171 strains of fungal
2. Claims 2-5,7-9 (partially) where they relate to probe 1415; a specific fungal probe for a fungal subject which includes yeasts but not *Candida* yeasts
3. Claims 1-5,7-9 (partially) where they relate to probe 1416; genus specific probe for *Candida* yeasts except. *C. lipolytica*
4. Claims 2-5,7-9 (partially) where they relate to probes 1G707 and 1859; specific for *C. lipolytica*
5. Claims 2-5,7-9 (partially) where they relate to probes 1814, 1816, 1813, 1860, 1812, 1857; probes defined as broadly inclusive for fungi
6. Claims 2-5,7-9 (partially) where they relate to probe 1858; probe specific for Zygomycetes